


Docking Stations

- DS4, DS6, DS8,
- DS4E, DS6E, DS8E

(with Ethernet module to work with the MedimaNet central infusion monitoring system)

- Available models accommodate 4, 6 or 8 syringe and volumetric pumps.
- Designed to save valuable space around the patient's bed.
- Allow fast and easy pumps mounting with automatic connection to the mains, without using interconnecting cables.



Technical data:

Power	100 - 240 V AC - 15%, +10%, 20 VA +35 VA for each pump
Alarm connectors	24V, 1A, NO/NC
Classification	Protection Class I
Electrical safety	EN60601-1, EN60601-1-1, EN60601-1-2, MDD 93/42/EEC - IIb
Dimensions in mm (W x H x D)	DS4/DS4E - 210 x 500 x 135 DS6/DS6E - 210 x 740 x 135 DS8/DS8E - 210 x 980 x 135
Weight (kg)	DS4/DS4E - 3.6 kg DS6/DS6E - 4.7 kg DS8/DS8E - 5.7 kg
Fastening	Vertical columns of the diameter within 16-32 mm, intended for attachment of heavy medical equipment.

Medima Sp. z o.o.
Al. Jerozolimskie 200, 02-486 Warsaw
tel. +48 (22) 313 22 66; fax +48 (22) 313 22 69
e-mail: medima@medima.com.pl
www.medima.com.pl

DSA Docking Station for Ambulances



- Allows fast and safe pump installation in ambulances (no additional cables needed).
- Enables automatic connection to DC power supply after installing.
- Protects the pump against vibration and impacts which may occur in an ambulance.

Technical data

Power supply	12–16 V DC, max 1 A with pump
Complies	EN60601-1, EN1789, MDD93/42/EEC-I
Mounting	To ambulance walls
Dimensions in mm (W x H x D)	372 x 180 x 158.5
Weight (kg)	1.6 (excl. connecting cords)

TP2 Duo Holder



- Allows fast connection and easy transport of two infusion pumps as one set
- Connection to the power supply with one cable
- Simple and light construction allows fast fastening to the IV pole

Technical data

Power supply	100 – 240, 50/60Hz, max 2 A
Dimensions in mm (W x H x D)	210 x 215 x 65
Weight (kg)	0.5

ST01 Medima Infusion Stand

For Infusion Pumps and Docking Stations

- Reliable structure ensures comfortable and safe fixing for docking stations and pumps.
- Enables easy and safe transport of the whole infusion set (docking station with a set of pumps - up to 8 pieces).

Technical data

Dimensions in mm (W x H x D)	480 x 1560 x 530
Weight in kg	15
Maximal load (kg)	35

Standard Infusion Stand

- Light and reliable structure for pumps mounting (up to three pieces) and transport.
- Available versions with or without power supply connection
- Regulated height within 130 – 225cm

Technical data

Dimensions in mm (W x H x D)	600 x 1950 x 600
Weight in kg	7
Maximal load (kg)	15

