

Safe and precise infusion

MEDIMA Line IV Sets




medima

MEDIMA Line IV Sets for Use with MEDIMA Volumetric Pumps




MEDIMA Line S

- Standard IV set for infusion of drugs and liquids
- Drop chamber, 20 drops/ml, filter 15µm
- Free flow protection clamp
- Roller clamp
- Luer lock
- Infusion port
- Length: 270 cm
- DEHP free
- Latex free



MEDIMA Line B

- IV set for infusion of blood and blood similar liquids
- Drop chamber, 20 drops/ml, filter 200µm
- Free flow protection clamp
- Roller clamp
- Luer lock
- Infusion port
- Length: 270 cm
- DEHP free
- Latex free



MEDIMA Line L

- IV set for infusion of light sensitive drugs, amber
- Drop chamber, 20 drops/ml, filter 15µm
- Free flow protection clamp
- Roller clamp
- Luer lock
- Infusion port
- Length: 270 cm
- DEHP free
- Latex free

MEDIMA Line IV sets can be used for administration of cytostatics

medima

Medima Sp. z o.o.

Al. Jerozolimskie 200, 02-486 Warsaw

tel. +48 (22) 313 22 66; fax +48 (22) 313 22 69

e-mail: medima@medima.com.pl

www.medima.com.pl